


Speed, Performance, Value - The Right Balance

- Large backlit LCD display
- Level indicator with adjustable levelling feet
- Colour coded keys
- Grade 304 stainless steel pan
- Simple user-friendly operation
- External calibration
- Zero Tracking
- AC Adapter
- Weighing
- Parts counting
- Percentage weighing
- Weight accumulation



Technical Specifications

Model	SCB 100	SCB 200	SCB 300	SCB 400	SCB 500	SCB 800	SCB 1200	SCB 2000	SCB 4000	SCB 6000	
Maximum capacity	100g	200g	300g	400g	500g	800g	1200g	2000g	4000g	6000g	
Readability	0.001g	0.01g	0.1g	0.01g	0.1g	0.1g	0.1g	0.1g	0.1g	1g	
Repeatability (S.D.)	0.003g	0.02g	0.2g	0.02g	0.2g	0.02g	0.2g	0.2g	0.2g	2g	
Draught shield dimensions	-							-			
Pan size	120mm & Oslash;							145mm & Oslash;			

General Specifications

Typical stabilisation time	2 Sec
Power supply	12VAC 50/60Hz 150mA adapter / 6 x AA batteries
Overall dimensions (wxdxh)	171 x 250 x 80mm
Net weight	0.9kg

Key Accessories

- ATP Thermal Printer
- Adam DU - Data Collection Program
- RS-232 to PC Cable



Adam Equipment - offices worldwide

www.adamequipment.com