

Wisconsin Oven Distributors, LLC.

W355 S9075 Godfrey Ln.,

Eagly, WI 53119

Phone: 262.594.3941 • Fax: 262.594.3939

E-Mail: tslaboven@centurytel.net • Web site: www.memmertbenchlabovens.com

Model IN55, Memmert Standard Incubators

[SPECIFICATIONS](#) · [STAINLESS STEEL INTERIOR](#) · [TEXTURED STAINLESS STEEL EXTERIOR](#) · [FURTHER DATA](#) · [PACKAGING DATA](#)

SPECIFICATIONS

Previous Models	INB400
NOTE	See downloadable PDF on right for detailed information on equipment, options and accessories.

STAINLESS STEEL INTERIOR

Volume	Approx. 1.9 ft ³ 53 L
Width (A)	15.7 in 400 mm
Height (B)	15.7 in 400 mm
Depth (Less 39 mm for Fan) (C)	13.0 in 330 mm
Number of Stainless Steel Grids (Standard Equipment)	1
Maximum Number of Grids	4
Maximum Loading per Grid	30 kg 66 lb
Maximum Loading of Chamber	80 kg 176 lb

TEXTURED STAINLESS STEEL EXTERIOR

Width (D)	23.0 in 585 mm
Height (Size 450, 750 with Castors) (E)	31.0 in 787 mm

Depth (with Door Handle), Door Handle + 56 mm (F)	20.2 in 514 mm
---	-------------------

FURTHER DATA

Voltage	230 V 115 V
Frequency	50/60 Hz
Power Consumption	1000 W
Working Temperature Range	Min. 5 °C (IN/INplus) 10 C (IF/IFplus) above ambient temperature up to +80 °C
Setting Temperature Range	+20 to +80 °C +68 to +176 °F
Setting Accuracy	0.1 °K

PACKAGING DATA

Approximate Net Weight	121 lb 55 kg
Approximate Gross Weight (packed in carton)	148 lb 67 kg
Approximate Height	37 in 94 cm
Approximate Width	28 in 70 cm
Approximate Depth	29 in 73 cm