

Titration

Application
support
available!

HI 903

Karl Fischer Titration System

HI 904

Coulometric Karl Fischer Titration System

Perfect for testing any parameter, HANNA's line of titrators combine an ultra-high precision titrant delivery system with optically regulated magnetic stirring, sophisticated endpoint determination, dynamic dosing and background drift correction algorithms.



Built-in magnetic stirrer



Field upgradable software via USB



- **HI 904 Coulometric Karl Fischer Titration System**
measures 000.5% to 100% water content with
precise dosing system capable of delivering as little as
0.125 µL of titrant with a single dose.

.....

- **HI 903 Karl Fischer Titration System**
measures 100% to 100 ppm water content.
Combines an ultra-high precision titrant delivery
system with optically regulated magnetic stirring,
sophisticated endpoint determination, dynamic dosing
and background drift correction algorithm

Perfect for

Food and Beverage • Pharmaceuticals • Cosmetics
Chemical and Petrol Manufacturing • Solvents