

MT-200ATC Protein Refractometer



The MT-200 is designed to measure the protein levels in blood serum. It works great for checking the protein levels in lobster blood to assure the quality is good.

Features:

- Three scales: serum protein, specific gravity, & refractive index
- Built-in automatic temperature compensation (ATC)
- Large easy to read lens

Specifications:

- Range:
 - Serum protein: 0-12g/100m
 - Specific gravity: 1.000-1.050sg
 - Refractive index: 1.330RI-1.360RI
- Resolution:
 - Serum Protein: 0.2g/100ml
 - Specific gravity: 0.005sg
 - Refractive index: 0.0005RI
- Weight: 0.6 lb