



Basic Info.

Model NO.:NSZ-806

Magnification:50-100X

Type:Stereomicroscope

Number of Cylinder:Monocular

Mobility:Desktop

Shape:Laminated Lenses

Principle:Optics

Illumination:Halogen Lamp

Export Markets:Global

Product Description

NSZ-806 Digital Zoom Stereo Microscope

Description:

Excellent Image Quality & Optical Performance with Infinity Parallel Optical System.

Comfortable Operation with Ergonomic Design Principle.

With LED light for both incident and transmitted illumination, Providing even illumination and life expectancy can reach 6000 hours.

A Wide Range of Accessories to be Equipped with the Microscope for Multi-purposelt is an Outstanding Assistant in

Industrial, Medical and Scientific Research Areas.

Specification:

Drawing Tube: Binocular Head, 20C Inclination

Eyepiece: EW10× /22

Zoom Objective: 0.8× -5×

0.8× -6.4×

0.8× -8×

Objective: Plan Achromatic Objective 1×

Zoom Ratio: 1: 6, 1: 8, 1: 10

Working Distance: 78mm

Focusing: Focusing Range: 105mm

Illumination: Transmission or Reflection LED Illumination, Brightness Adjustable