EBA 200 TABLETOP CENTRIFUGE

PRESS 'START' IN CASE OF EMERGENCY



Compact clinical tabletop centrifuge

High-performance, built-in, 8-place rotor with compact frame and all-metal containment for added safety. Ideal for blood, urine and pediatric tubes.

STAT, PPP & Coagulation

Obtain consistent platelet counts less than 10,000 platelets/µL in 3 minutes or less (model EBA 200S)

Run-Silent™ Operation

Excessive noise in the laboratory can lead to distractions, stress, anxiety, high blood pressure and chronic fatigue. Our centrifuges are as quiet as a normal conversation and will not increase the noise levels within your lab.

Easy-Lift™

Our Easy-Lift Lid, featuring stay-in-place positioning and lid drop protection, is ergonomically designed to limit repetitive stress injuries, spillage and accidental sample disruption when loading and unloading your centrifuge.

Safety-Lock™

Our Safety-Lock Lids are self-closing, auto-locking and come with either single and dual steel lid locks. This ensures a safe, tight closure of the centrifuge chamber during operation. The centrifuge will not operate if lid is not secure. For added safety, our centrifuges feature all-metal containment and adhere to all national and international safety standards and certifications.

Five (5) Year Warranty

FDA registered centrifuge manufacturer









HETTICH EBA 200 CENTRIFUGE

APPLICATION PACKAGES

All Hettich centrifuge packages come complete with everything you need to perform your application – centrifuge, rotor, buckets, adapters, inserts, bio-containment lids (optional where applicable), five (5) year warranty and free support.

EBA 200 BLOOD TUBE PACKAGE - 8

and the second	
-	

BE PACKA	GE - 8	EBA 200	EBA 200S	EBA 200	EBA 200S	
Tube Qty	Volume (mL)	Size (mm)	R	PM	G-F	ORCE
8	1.6-5	13 x 75	6,000	8,000	2,697	4,794
8	4-7	13 x 100	6,000	8,000	3,461	6,153
8	8.5-10	16 x 100	6,000	8,000	3,461	6,153
4	8	16 x 125	6,000	8,000	3,461	6,153

PKG. NO. EBA 200	
200BL00D8	
PKG. NO. EBA 200S	
200SBL00D8	

EBA 200 PEDIATRIC TUBE PACKAGE - 8



Tube Qty Volume (mL) Size (mm) RPM G-FORCE 8 0.5 10.7 x 36 6.000 8.000 2.214 3.93	TUBE PACKAGE - 8		EBA 200	EBA 200S	EBA 200	EBA 200S	
8 0.5 10.7 x 36 6.000 8.000 2.214 3.93	Tube Qty	Volume (mL) Size (mm) RPM		G-FORCE			
7	8	0.5	10.7 x 36	6,000	8,000	2,214	3,935

PKG. NO. EBA 200
200BL00D8-PED
PKG. NO. EBA 200S
200SBL00D8-PED

EBA 200 URINALYSIS PACKAGE - 4



Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE
4	12 (KOVA)	17 x 102	6,000	3,461

EBA 200 CUSTOM PACKAGE



Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PKG. NO.

Additional packages and configurations are available for this device, please inquire.

We offer a complete centrifuge product range from small microliter centrifuges, compact tabletop centrifuges, medium to large benchtop centrifuges, high capacity benchtop and floor standing centrifuges for general purpose and specialty applications, inquire for more details.

→ Learn more at: www.hettweb.com/eba200

DISTRIBUTED BY:

