

**【RE100-Pro  
LCD Digital Rotary Evaporator】**

**Features**

- Easy to read large digital LCD screen displays heating temperature, rotation speed and timing
- Speed range from 20 to 280 rpm
- Large 5L water-oil heating bath with heating temperature range of RT to 180°C
- Heating bath with precise temperature control and adjustable safety circuit
- Designed a detachable operating panel allows remote control
- Patented condenser(cooling surface 1,500 cm<sup>2</sup>) with excellent cooling effect
- Motorized lift with quick-action, automatic release evaporating flask to top position in case

