

Rotary Microtome

Rotary Microtome YD-1508A



Product Description

Paraffin Microtome YD-1508A(B), from China professional Microtome manufacturer Jinhua YIDI Medical Appliance Co., Ltd

The machines have the characteristics of the good stability, high precision, operating easily, working continually and evenly.

Technical Data:

- 1) Slice thickness range: 1-30 micron(YD -1508A), 1-25micron(YD-1508B)
- 2) Least slice adjusting graduation: 1micron
- 3) Slice precision: $\pm 5\%$
- 4) Maximal slice section: 50× 30mm (40× 30mm)