

QTS-Smart series



VT-13QTS Ultrasonic Cleaner 13L

Products Introduction:

Features:

- Turbo power: Normal and soft function to resolve the cleaning blind spot .
 - Auto-degas function for laboratory purposes
 - Memory function
 - Automatically ultrasonic frequency tracking
 - LED display for temperature and timer
 - Timer:1-99minutes
 - Temperature:0-80°C
 - POSCO SUS304 tank
 - User-friendly and clear panel
-
- High-performance transducer systems

Ultrasonic cleaning:

The ultrasonic cleaning process provides high quality cleaning and exceeds, by ease of use and end result, conventional cleaning methods. The tiny microscopic bubbles generated by the transducers provide cleaning of very difficult, irregular and complex parts. Even inaccessible areas can be cleaned without any trace of physical cleaning (scratches, etc...).

Ultrasonic cleaner is widely applied to different fields:dental instruments, lab instruments,hospital,instruments,

blind,PC board,hardware,chemical industry,circuit board,carburetor,printer head,and etc.

Specifications:

Model	VT-13QTS
Voltage	AC100~120V, 60Hz AC220~240V, 50Hz
Frequency	40KHz
Power	300W
Tank capacity	13L
Tank size(mm)	330×300×150 (L×W×H)
Overall size(mm)	360×330×310 (L×W×H)
Unit weight	11.8KG
Lid	Have
Heating	Have
Drainage	Have
Cleaning basket	Optional