

## QTS-Smart series



### VT-30QTS ultrasonic cleaner 3L

#### Products Introduction:

#### Features:

- \*Turbo power: Normal and soft function to resolve the cleaning blind spot .
- \*Degas function for laboratory purposes
- \* Memory function
- \*Automatically ultrasonic frequency tracking
- \* LED display for temperature and timer
- \*Timer: 1-99minutes
- \*Temperature:0-80°C
- \*POSCO SUS304 tank
- \*User-friendly and clear panel
- \*High-performance transducer systems

#### Ultrasonic cleaning:

The ultrasonic cleaning process provides high quality cleaning and exceeds, by ease of use and end result, conventional cleaning methods. The tiny microscopic

bubbles generated by the transducers provide cleaning of very difficult, irregular and complex parts. Even inaccessible areas can be cleaned without any trace of physical cleaning (scratches, etc...).

Ultrasonic cleaner is widely applied to different fields:dental instruments, lab instruments,hospital,instruments,

blind,PC board,hardware,chemical industry,circuit board,carburetor,printer head,and etc.

### **Specifications:**

Model	VT-30QTS
Voltage	AC100~120V,60Hz AC220~240V,50Hz
Frequency	40KHz
Power	100W
Tank capacity	3L
Tank size(mm)	240×140×100 (L×W×H)
Overall size(mm)	270×170×240 (L×W×H)
Unit weight	3.8KG
Lid	Have
Heating	Have
Drainage	None
Cleaning basket	Optional