

# Ultrasonic Sonicator



## VT650-IIN, VT900-IIDN Ultrasonic Homogenizer

Model VT650-IIN, VT900-IIDN ultrasonic homogenizer are multi-function and multipurpose instruments that utilize strong super sound produce empty effect of taking among liquid, the instruments of multi-function, multiuse of dealing with the super sound to the material, can be used in the brokenness of many kinds of animal and plant cells , virus cell, at the same time , can use emulsification , separates , melts , draws , subdues steeping , washes and accelerates the chemical reaction etc. evenly. Widely used in such fields as biochemistry, microbiology, medicine chemistry, surface chemistry, physics, zoology, etc.

### Main Technical Parameters:

MODEL	VT650-IIN	VT900-IIDN
Frequency	20-25 KHz	20-25 KHz
Power	650 W	900 W
Scessories	6mm	6mm
Optional Tip	2,3,10,12,15mm	2,3,10,12,15mm
Crushing Capacity	0.5-500 ml	0.5-500 ml
Duty Ratio (%)	1-99 %	0.1-99.9%
Supply	220/110V 50Hz/60Hz	220/110V 50Hz/60Hz
Power supply chassis Size and Weight	40*28*22(cm);12.1Kg	40*28*22(cm);12.1Kg
Mainframe Weight	14 Kg	14.2 Kg
Packing Size	534*295*435mm	534*295*435mm