



Product Description

Features:

Highly Bright Optical System

Mechanical Tube Length 160mm

Monocular Head Inclined at 45°, 360° Rotatable

Compensation Free Binocular Head Inclined at 30°, 360° Rotatable, Interpupillary Distance 55-57mm

Compensation Free Trinocular Head Inclined at 30°, 360° Rotatable, Interpupillary Distance 55-57mm

Sliding Binocular Head Inclined at 45°, 360° Rotatable, Interpupillary Distance 55-57mm

Double Layers Mechanical Stage:

Size 140x160mm, Moving range 75x50mm

Coaxial Coarse and Fine Focus Adjustment Focusing Range 30mm, division of FA 0.002mm

Coarse Adjustment Tension Control with Focusing Stop to Prevent Objectives and Slide from Damaging

Rack and Pinion Adjustment for condenser

Built-in Illumination with Halogen Lamp, Adjustable Brightness

A Wide Range of Accessories to be Equipped with this Microscope for Multi-purpose