



CENTRIFUGES

EBA 270

TABLETOP CENTRIFUGE

90° SST, easy as 1, 2, 3.

HORIZONTAL CLINICAL TABLETOP CENTRIFUGE

90° Built-in, swing-out rotor for flat separation of all standard blood, gel, urine and pediatric tubes. Ideal for physicians and veterinary offices.

RUN-SILENT™ OPERATION

Excessive noise in the laboratory can lead to distractions, stress, anxiety, high blood pressure and chronic fatigue. Our centrifuges are as quiet as a normal conversation and will not increase the noise levels within your lab.

EASY-LIFT™

Our Easy-Lift Lid, featuring stay-in-place positioning and lid drop protection, is ergonomically designed to limit repetitive stress injuries, spillage and accidental sample disruption when loading and unloading your centrifuge.

SAFETY-LOCK™

Our Safety-Lock Lids are self-closing, auto-locking and come with either single or dual steel lid locks. This ensures a safe, tight closure of the centrifuge chamber during operation. The centrifuge will not operate if lid is not secure. For added safety, our centrifuges feature all-metal containment and adhere to all national and international safety standards and certifications.

FIVE (5) YEAR WARRANTY

FDA REGISTERED MANUFACTURER



hettweb.com/eba270



MAX CAPACITY
6 x 15 mL



MAX RPM/RCF
4,000/2,254



DIMS. (HxWxD)
9.5" x 12.9" x 15.3"




TUBE TYPES
BLOOD, PEDIATRIC, URINE

HETTICH EBA 270 APPLICATION PACKAGES




EBA 270 **COMPLETE** PACKAGE - 6


	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PKG. NO.
	6	1.6 – 7	13/16 x 75	4,000	1,807	270COMPLETE
	6	4 – 10	13/16 x 100	4,000	2,254	
	6	0.5 (Pediatric)	10.7 x 36	4,000	1,359	
	6	12 (KOVA)	17 x 102	4,000	2,254	

All Hettich centrifuge packages come complete with everything you need to perform your application – centrifuge, rotor, buckets, adapters, inserts, bio-containment lids (optional where applicable), five (5) year warranty & free support.





SMALL IN SIZE, BIG ON RESULTS
Blood separation on the EBA 270 spinning a standard 15 mL Gel Tube.





hettweb.com/eba270

CONTACT YOUR LOCAL HETTICH REPRESENTATIVE TO LEARN MORE:



Discover more solutions online at: hettweb.com

HETTICH LAB TECHNOLOGY | North America **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com

