

Motorized Pipette Filler

Levo Plus



Features:

- ❑ A full volume range of 1-100ml
- ❑ Easy to one-handed operation
- ❑ Lightweight, ergonomic design allows longer, fatigue-free pipetting
- ❑ Suitable for all types of pipettes 1-100ml
- ❑ Lithium-ion battery provides long runtime charge and no battery memory problems
- ❑ Low battery indication
- ❑ Powerful motor fills a 25ml pipette in <5 seconds
- ❑ Supplied with integral but replaceable 0.45um filter
- ❑ Can be used while recharging

Main Data Description	
Aspirate Speeds	8
Dispense Speeds	8 & Gravity Dispense
Battery	Lithium-Ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	Approx.2 Hours
Pipette Types	Glass or Plastic Pipettes(0.1-100ml); Pasteur Pipettes
Filter	0.45µm Hydrophobic

* O.E.M. option also available.
* Design & Specification of Pipette Filler & its accessories are subject to change without prior notice.

Ordering information

Cat No.	Description
740200010000	Levo Plus with AC adapter, spare 0.45um filter and wall stand, Euro plug
740200020000	Levo Plus with AC adapter, spare 0.45um filter and wall stand, North America plug

Accessories

Cat No.	Description
17000103	Filter 0.45um, 1pc each pk

Dragon Laboratory Instruments Limited
F3, Block 12, No.28 Yuhua Road, Shunyi District, Zone B,
Tian Zhu Airport Industrial Park, Beijing, 101318 P. R. China
Tel: 86-10-85653451 85653452 85653453
Fax: 86-10-85653383
E-mail: info@dragon-lab.com
Service: service@dragon-lab.com
Service Hotline: 400-1100-556
Website: [Http://www.dragon-lab.com](http://www.dragon-lab.com)

US office
DRAGONLAB, LLC
500 Four Rod Road, Suite 122, Berlin, CT 06037 USA
Tel: 1-860-828-5289
Fax: 1-860-828-5389
E-mail: info@dragonlabusa.com
Website: [Http://www.dragon-lab.com](http://www.dragon-lab.com)