

【MS-H280-Pro  
LED Digital Magnetic Hotplate Stirrer】

**Features of MS-H280-Pro**

**Digital control system is ideal for applications for physical - chemical analysis and biotech labs.**

- Heating temperature up to 280°C
- Stainless steel with ceramic coated hotplate resists to corrosion
- External temperature sensor PT1000 is available and control accuracy of  $\pm 1^{\circ}\text{C}$ , real-time control medium temperature.
- LED display temperature and speed.
- 50°C HOT warning indicates residual hotplate temperature.
- ABS casing with fire retardant material, resistance of weak acid and alkali



MS-H280-Pro  
with external temperature sensor PT1000



MS-H280-Pro  
stainless steel with ceramic coated hotplate

