

**【OS40/20-Pro
LCD Digital Overhead Stirrer】**

Overhead stirrer is suitable for mixture of high viscosity liquid or solid-liquid. Widely used in chemical synthesis, pharmaceutical, physical and chemical analysis, petrochemical, cosmetics, health care, food, biotechnology and other fields.



OS40-Pro
LCD Digital Enhanced Overhead Stirrer
OS20-Pro
LCD Digital Overhead Stirrer



OS40-Pro/ OS20-Pro
with universal plate stand

Features of OS40-Pro/ OS20-Pro

- High reliability, long continuous service, outstanding performance
- Cover the speed range from 50 to 2200rpm, the speed is accurate and controlled
- Brushless DC motor drive, low noise and maintenance free
- Setting value and actual value can be displayed, real-time monitoring torque and speed via LCD display
- Overload protection, motor would automatically stop in case of continuous overload
- Anti-spill samples, PC control system starts up smoothly, prevent samples overflow
- Extended application, PC can control and document all operation values via RS232 interface.

