

**【SK-R330-Pro
LCD Digital Rocking Shaker】**

Rocking shaker is suitable for mixing applications in wide fields including, gel staining and destaining, washes, hybridisations, immune precipitation and Western blots etc.



**SK-R330-Pro
LCD Digital Rocking Shaker
Attachment with non-slip mat**

Features of SK-R330-Pro

- Smooth see-saw rocking motion with fixed angle 9°
- Load up to 10kg
- Speed range of 10-70rpm
- DC brushless motor provides reliable and quiet running, maintenance free
- Designed for continuous or timed operation
- Digital speed control and LCD display, speed and timer can be viewed simultaneously
- Wide range of accessory platform to be used for variable size vessels, and non-slip mat holds vessels prevent slipping

